

**Note**


Unless otherwise indicated, the left-hand side of the image faces the middle ear.



**ARMSTRONG BEVELED GROMMETS**

Material	Cat. No.	ID	Description
Fluoroplastic	14-0242	1.14 mm	
Fluoroplastic	*14-5700	1.14 mm	
Fluoroplastic	7014-5061	1.14 mm	Blue
Fluoroplastic	*7014-5760	1.14 mm	Blue
Fluoro-PacifiC	7014-0267	1.14 mm	
Fluoro-PacifiC	**7014-3067	1.14 mm	
Fluoro-PacifiC	*7014-6067	1.14 mm	
Fluoro-PacifiC	7014-0395	1.14 mm	Arm. V Style, Blue
Fluoro-PacifiC	**7014-0396	1.14 mm	Arm. V Style, Blue
Fluoro-PacifiC	*7014-0397	1.14 mm	Arm. V Style, Blue
Silicone	24-0050	1.14 mm	w/Tab
Ultrasil	7024-1050	1.14 mm	w/Tab
Ultrasil	**7024-1350	1.14 mm	w/Tab
Ultrasil	*7024-5755	1.14 mm	w/Tab

Designed by B.W. Armstrong, M.D., Charlotte, NC.



**ARMSTRONG BEVELED PLAIN END VENTILATION TUBES**

Material	Cat. No.	ID	Length
Fluoroplastic	14-0244	1.14 mm	7.5 mm
Fluoroplastic	*14-5706	1.14 mm	7.5 mm
Fluoro-PacifiC	7014-0345	1.14 mm	7.5 mm
Fluoro-PacifiC	**7014-0346	1.14 mm	7.5 mm
Fluoro-PacifiC	*7014-0347	1.14 mm	7.5 mm
Fluoro-PacifiC	7014-0398	1.14 mm	Arm V, 4.5 mm
Fluoro-PacifiC	*7014-0399	1.14 mm	Arm V, 4.5 mm
Fluoro-PacifiC	**7014-0400	1.14 mm	Arm V, 4.5 mm
Fluoro-PacifiC	7014-0401	1.14 mm	Arm V, 7.5 mm
Fluoro-PacifiC	*7014-0402	1.14 mm	Arm V, 7.5 mm
Fluoro-PacifiC	**7014-0403	1.14 mm	Arm V, 7.5 mm
Silicone	24-0052	1.14 mm	7.5 mm
Ultrasil	7024-1052	1.14 mm	7.5 mm

Designed by B.W. Armstrong, M.D., Charlotte, NC.



**BAXTER BUTTON VENTILATION TUBE**


Material	Cat. No.	ID	Description
Fluoro-PacifiC	7014-0271	97 mm	Blue
Fluoro-PacifiC	**7014-0273	97 mm	Blue
Fluoro-PacifiC	*7014-0275	97 mm	Blue
Fluoro-PacifiC	7014-0270	1.27 mm	Blue
Fluoro-PacifiC	**7014-0272	1.27 mm	Blue
Fluoro-PacifiC	*7014-0274	1.27 mm	Blue



**BERGER "V" BOBBIN VENTILATION TUBE**

Material	Cat. No.	ID	Description
Titanium	14-5282	1.27 mm	

Designed by Eric S. Berger, M.D., Aberdeen, WA. U.S. Patent No. 4,080,171.



**BUTTERFLY VENTILATION TUBE**

Material	Cat. No.	ID	Description
Silicone	24-0073	1.27 mm	5.3mm L, 7.6mm W
Silicone	**7024-5018	1.27 mm	5.3mm L, 7.6mm W
Silicone, Ultrasil	7024-1073	1.27 mm	5.3mm L, 7.6mm W



**CANALE VENTILATION TUBE**

Material	Cat. No.	ID	Description
Silicone	24-0080	1.14 mm	Blue

Developed in conjunction with James Canale, M.D., Blytheville, AR.



**COLLAR BOBBIN® VENT TUBES**

Material	Cat. No.	ID	Description
Titanium	14-5287	1.14 mm	
Fluoroplastic	14-5288	.75 mm	
Fluoroplastic	14-5289	1.14 mm	
Titanium	14-5297	1.27 mm	
Fluoroplastic	14-5298	1.27 mm	

Designed for Bruce Black, M.D., Brisbane, Australia.



**COLLAR BUTTON VENTILATION TUBES**


Material	Cat. No.	ID	Description
Fluoroplastic	14-5214	1.27 mm	White
Fluoroplastic	**14-5004	1.27 mm	White
Fluoroplastic	*14-5705	1.27 mm	White
Fluoroplastic	*7014-5062	1.27 mm	Blue
Fluoroplastic	*7014-5060	1.27 mm	Blue
Fluoroplastic	14-5217	1.27 mm	White w/Wire
Fluoro-PacifiC	7014-5314	1.27 mm	White
Fluoro-PacifiC	**7014-5316	1.27 mm	White
Fluoro-PacifiC	*7014-5317	1.27 mm	Blue
Silicone	24-0014	1.27 mm	Blue
Silicone, Ultrasil	7024-1014	1.27 mm	Blue
Silicone, Ultrasil	**7024-1314	1.27 mm	Blue
Silicone, Ultrasil	*7024-1060	1.27 mm	Blue



**DONALDSON VENTILATION TUBES**

Material	Cat. No.	ID	Description
Silicone	24-0015	1.14 mm	
Soft Silicone	24-0060	1.14 mm	
Silicone	*24-5702	1.14 mm	
Silicone	24-0016	1.14 mm	w/Tab
Silicone	**7024-5012	1.14 mm	w/Tab
Silicone	*24-5707	1.14 mm	w/Tab
Fluoroplastic	14-6043	1.14 mm	
Ultrasil	7024-1015	1.14 mm	

Designed for James A. Donaldson, M.D., Seattle, WA.



**DURAVENT™ TUBES**

Material	Cat. No.	ID	Description
Silicone	24-0075	1.27 mm	1.37 mm IF, Blue
Silicone	*24-5775	1.27 mm	1.37 mm IF, Blue
Ultrasil	7024-1075	1.27 mm	1.37 mm IF, Blue
Ultrasil	**7024-1375	1.27 mm	1.37 mm IF, Blue

Designed by Bruce Black, M.D., Brisbane, Australia.



**DUTCHER VENTILATION TUBE**

Material	Cat. No.	ID	Description
Silicone, Ultrasil	7024-0054	1.27 mm	4.32 mm Medial Flange Width

Designed for Paul Dutcher, M.D., Rochester, NY.



**FEUERSTEIN SPLIT TUBE**

Material	Cat. No.	ID	Length
Fluoroplastic	14-0248	1.00 mm	9 mm

Designed for Sidney S. Feuerstein, M.D., New York, NY.



**JAHN HYDROXYLVENT® TUBES**

Material	Cat. No.	ID	Description
Hydroxylapatite	14-0906	1.00 mm	9.5 mm Long
Hydroxylapatite	14-0926	1.27 mm	9.5 mm Long
Hydroxylapatite	14-0927	1.00 mm	w/Fpl. Liner

Packaged individually. Designed for Anthony Jahn, MD, New York, NY.



**LINDEMAN-SILVERSTEIN ARROW VENTILATION TUBE**

Material	Cat. No.	ID	Description
Fluoroplastic	14-0240	1.14 mm	

Designed for Roger Lindeman, M.D., Seattle, WA, and Herbert Silverstein, M.D., Sarasota, FL.



**MICRON® TITANIUM BOBBIN VENTILATION TUBES**


Material	Cat. No.	ID	Description
Titanium/Fluoro	14-5220	1.00 mm	w/Flange Holes
Titanium/Fluoro	14-5271	1.27 mm	
Titanium/Fluoro	14-5272	1.52 mm	
Titanium	14-5221	0.76 mm	Infant
Titanium	14-5215	1.00 mm	w/Flange Holes
Titanium	14-5219	1.00 mm	"Unibobbin" (FH on 1 side)
Titanium	14-5280	1.00 mm	
Titanium	14-5279	1.14 mm	
Titanium	14-5273	1.27 mm	w/Spoon Flange
Titanium	14-5276	1.27 mm	+ Hydromer™ Coat
Titanium	*7014-5776	1.27 mm	+ Hydromer™ Coat
Titanium	*14-5781	1.27 mm	
Titanium	14-5281	1.27 mm	
Blue Titanium	7014-5381	1.27 mm	Blue
Blue Titanium	7014-5881	1.27 mm	Blue, 12/box



**MINI-REVERSIBLE VENTILATION TUBE**

Material	Cat. No.	ID	Description
Silicone	24-0042	0.89 mm	Reversible

Designed for Burton J. Cohen, M.D., Louisville, KY.



**MODIFIED T-TUBES**

Material	Cat. No.	ID	Length
Silicone	24-0072	1.32 mm	4.80 mm
Silicone	**7024-5017	1.32 mm	4.80 mm
Silicone	*24-5704	1.32 mm	4.80 mm
Soft Silicone	24-0074	1.32 mm	4.80 mm
Ultrasil	7024-1072	1.32 mm	4.80 mm



**MORETZ TYPE VENTILATION TUBES**

Material	Cat. No.	ID	Description
Titanium	14-5265	1.27 mm	
Polyethylene	7014-5498	1.27 mm	Blue
Polyethylene	**7014-3498	1.27 mm	Blue
Polyethylene	*7014-5798	1.27 mm	Blue



**PER-LEE VENTILATION TUBE**

Material	Cat. No.	ID	Description
Silicone	24-0026	1.57 mm	60°, 15 mm Long

Designed for John H. Per-Lee, M.D., Emory Univ., Atlanta, GA.



### PAPARELLA-TYPE VENTILATION TUBES

Material	Cat. No.	ID	Description
Fluoroplastic	14-0252	1.00 mm	w/Tab
Fluoroplastic	14-0254	1.27 mm	w/Tab
Silicone	24-0043	1.00 mm	w/o Tab
Silicone	24-0044	1.00 mm	w/Tab
Silicone	**7024-5016	1.00 mm	w/Tab
Silicone	**24-5701	1.00 mm	w/Tab
Silicone	7014-0280	1.14 mm	w/Tab
Silicone	24-0046	1.27 mm	w/Tab
Silicone	**7024-3046	1.27 mm	w/Tab
Silicone	**7024-5746	1.27 mm	w/Tab
Silicone	24-0048	2.00 mm	w/Tab
Soft Silicone	24-0058	1.00 mm	w/Tab
Soft Sil-Ultrasil	7014-0281	1.14 mm	w/Tab
Soft Sil-Ultrasil	**7024-1380	1.14 mm	w/Tab
Soft Sil-Ultrasil	**7024-5780	1.14 mm	w/Tab
Silicone, Ultrasil	7024-1044	1.00 mm	w/Tab
Silicone, Ultrasil	**7024-1344	1.00 mm	w/Tab
Silicone, Ultrasil	7024-0280	1.14 mm	w/Tab
Silicone, Ultrasil	**7024-0281	1.14 mm	w/Tab
Silicone, Ultrasil	**7024-0282	1.14 mm	w/Tab
Silicone, Ultrasil	7024-1046	1.27 mm	w/Tab

Designed for Michael Paparella, M.D., Minneapolis, MN.



### POPE BEVELED GROMMET VENTILATION TUBES

Material	Cat. No.	ID	Description
Polyethylene	14-5250	1.14 mm	Blue
Polyethylene	**7014-3250	1.14 mm	Blue
Polyethylene	**7014-5750	1.14 mm	Blue


Designed for T.H. Pope, Indianapolis, IN.



### REUTER BOBBIN VENTILATION TUBES

Material	Cat. No.	ID	Description
Fluoroplastic	14-5212	1.00 mm	1 mm IF, w/o Holes
Fluoroplastic	*14-5709	1.00 mm	1 mm IF, w/o Holes
Fluoroplastic	14-5213	1.14 mm	1 mm IF, w/o Holes
Fluoroplastic	**14-5019	1.14 mm	1 mm IF, w/o Holes
Fluoroplastic	14-5224	1.14 mm	1 mm IF, w/Wire
Fluoroplastic	14-5225	1.14 mm	1 mm IF
Fluoroplastic	*14-5710	1.14 mm	1 mm IF
Fluoroplastic	14-5227	1.27 mm	0.5 mm IF
Fluoro-PacifiC	7014-5319	1.00 mm	1 mm IF
Fluoro-PacifiC	**7014-5320	1.00 mm	1 mm IF
Fluoro-PacifiC	*7014-5321	1.00 mm	1 mm IF
Fluoro-PacifiC	7014-5326	1.00 mm	1 mm IF
Fluoro-PacifiC	**7014-5327	1.00 mm	1 mm IF
Fluoro-PacifiC	*7014-5328	1.00 mm	1 mm IF
Stainless Steel	14-5210	1.00 mm	1 mm IF
Stainless Steel	*14-5701	1.00 mm	1 mm IF
Stainless Steel	14-5211	1.00 mm	Satin, 1 mm IF W/Wire
Stainless Steel	14-5216	1.00 mm	Satin, 1 mm IF
Stainless Steel	14-5218	1.00 mm	Satin, 1 mm IF
Stainless Steel	14-5208	1.14 mm	1 mm IF
Stainless Steel	14-5209	1.27 mm	1 mm IF
Titanium	14-5215	1.00 mm	Micron Finish

Designed for S. Harold Reuter, M.D., Houston, TX.




### RICHARDS T-TUBES

Material	Cat. No.	ID	Length	Width
Silicone	24-0071	1.14 mm	12 mm	9.8 mm
Silicone	**7024-5014	1.14 mm	12 mm	9.8 mm
Soft Silicone	24-0059	1.14 mm	12 mm	9.8 mm
Silicone	*24-5706	1.14 mm	12 mm	9.8 mm
Ultrasil	7024-1071	1.14 mm	12 mm	9.8 mm



### ROCK PEDIATRIC VENTILATION TUBE

Material	Cat. No.	ID	Description
Fluoroplastic	14-0086	1.00 mm	0.89 mm IF



### RUBE VENTILATION TUBE

Material	Cat. No.	ID	Length	Width	Description
Silicone	7014-0260	1.14 mm	4.5 mm	6.6 mm	w/Tab



### SHAH VENTILATION TUBES

Material	Cat. No.	ID	Description
Fluoroplastic	14-2230	0.76 mm	w/o Wire
Fluoroplastic	**14-5020	0.76 mm	w/o Wire
Fluoroplastic	14-6044	1.14 mm	w/o Wire
Fluoroplastic	**14-5001	1.14 mm	w/o Wire
Fluoroplastic	14-0005	1.14 mm	w/Wire
Fluoroplastic	**14-5015	1.14 mm	w/Wire

Designed for Mr. N. Shah, FRCS, DIO, London, UK.



### SHEA PARASOL VENTILATION TUBES

Material	Cat. No.	ID	Description
Silicone	24-0035	1.00 mm	w/Tab, Small Flange
Silicone	24-0036	1.00 mm	w/Tab
Silicone	*24-5705	1.00 mm	w/Tab
Silicone	24-0038	1.50 mm	w/Tab
Ultrasil	7024-1036	1.00 mm	w/Tab

Designed for John J. Shea, M.D., Memphis, TN. U.S. Patent No. 3,871,380.



### SHEA VENTILATION TUBE

Material	Cat. No.	ID	Description
Silicone	24-0034	0.76 mm	w/Tab

Designed for John J. Shea, M.D., Memphis, TN. U.S. Patent No. D239,330.



### SHEEHY-TYPE COLLAR BUTTON


Material	Cat. No.	ID	Description
Fluoroplastic	7014-5971	1.27 mm	
Fluoroplastic	**7014-3071	1.27 mm	
Fluoroplastic	*7014-5771	1.27 mm	
Fluoro-PacifiC	7014-1080	1.27 mm	White
Fluoro-PacifiC	**7014-3075	1.27 mm	White
Fluoro-PacifiC	*7014-5780	1.27 mm	White



### SHEPARD GROMMET VENTILATION TUBES

Material	Cat. No.	ID	Description
Fluoroplastic	14-0007	1.02 mm	w/Tab
Fluoroplastic	14-0008	1.14 mm	w/Tab
Fluoroplastic	14-0027	1.14 mm	
Fluoroplastic	**14-5003	1.14 mm	
Fluoroplastic	14-0028	1.00 mm	w/Wire
Fluoroplastic	**14-5009	1.00 mm	w/Wire
Fluoroplastic	*14-5702	1.00 mm	w/Wire
Fluoroplastic	14-0029	1.00 mm	
Fluoroplastic	**14-5002	1.00 mm	
Fluoroplastic	*14-5711	1.00 mm	
Fluoroplastic	14-0030	1.14 mm	w/Wire
Fluoroplastic	**14-5010	1.14 mm	w/Tab and Wire
Fluoroplastic	*14-5704	1.14 mm	w/Tab and Wire
Fluoroplastic	7014-1030	1.02 mm	
Fluoroplastic	**7014-1031	1.02 mm	
Fluoroplastic	**7014-1032	1.02 mm	
Silicone	24-0020	1.00 mm	w/Tab
Silicone	**7024-5008	1.00 mm	w/Tab
Polyethylene	14-6028	1.00 mm	w/Tab and Wire, Blue
Polyethylene	14-6029	1.00 mm	w/Tab, Blue
Polyethylene	14-6030	1.14 mm	w/Tab and Wire, Blue
Polyethylene	14-6031	1.14 mm	w/Tab, Blue
Ultrasil	7024-1020	1.00 mm	


Designed for Marvin G. Shepard, M.D., Dallas, TX.



### SPOON BOBBIN VENTILATION TUBES

Material	Cat. No.	ID	Description
Fluoroplastic	*7014-5224	1.00 mm	w/Wire, w/o Holes
Stainless Steel	14-5240	1.00 mm	w/o Wire
Stainless Steel	14-5241	1.00 mm	w/Wire
Stainless Steel	**7014-5225	1.00 mm	w/Wire
Titanium	14-5273	1.27 mm	w/Wire, w/o Holes

Designed for Baldev K. Devgan, M.D., Malibu, CA.



### STRAIGHT VENTILATION TUBES


Material	Cat. No.	ID	Description
Fluoroplastic	14-0032	0.89 mm	12 mm Long
Fluoroplastic	**14-5010	0.89 mm	12 mm Long
Fluoroplastic	14-0033	1.14 mm	12 mm Long
Fluoroplastic	**14-5007	1.14 mm	12 mm Long
Fluoroplastic	*14-5703	1.14 mm	12 mm Long
Fluoroplastic	14-0034	0.89 mm	7 mm Long
Fluoroplastic	*14-5707	0.89 mm	7 mm Long
Fluoro-PacifiC	7014-1035	0.89 mm	7 mm Long
Fluoro-PacifiC	**7014-1036	0.89 mm	7 mm Long
Fluoro-PacifiC	*7014-1037	0.89 mm	7 mm Long
Fluoro-PacifiC	14-0035	1.14 mm	7 mm Long
Fluoro-PacifiC	**14-5011	1.14 mm	7 mm Long
Fluoro-PacifiC	*14-5708	1.14 mm	7 mm Long
Fluoro-PacifiC	7014-1039	1.14 mm	7 mm Long
Fluoro-PacifiC	**7014-1040	1.14 mm	7 mm Long
Fluoro-PacifiC	*7014-1041	1.14 mm	7 mm Long
Polyethylene	14-0016	0.86 mm	7 mm Long, Clear
Polyethylene	*14-5713	0.86 mm	7 mm Long, Clear
Polyethylene	14-6033	1.14 mm	12 mm Long, Blue
Polyethylene	14-6035	1.14 mm	7 mm Long, Blue



### SULTAN VENTILATION TUBES

Material	Cat. No.	ID	Description
Fluoroplastic	14-5230	1.27 mm	w/Tab
Fluoroplastic	**14-5005	1.27 mm	w/Tab

Designed for Andre A. Sultan, M.D., France.



### TINY T-TUBES

Material	Cat. No.	ID	Length	Width
Silicone	7024-0076	1.32 mm	4.0 mm	4.0 mm
Ultrasil	7024-1076	1.32 mm	4.0 mm	4.0 mm

Designed in conjunction with Richard L. Goode, MD, Palo Alto, CA.



### VENTURI+PLUS VENTILATION TUBE

Material	Cat. No.	ID	Description
Titanium	7014-5278	1.27 mm	Hydromer Coated

#### Notes

All ventilation tubes are packed one per package, six per box, unless otherwise specified.

\* Twinpack: two vent tubes per package, 30 pkgs./box, 60 vent tubes total

\*\* 30 Pack (Singles).

• Twelvepack.

#### Trademarks

DuraVent is a trademark of Gyrus ENT LLC.

Collar Bobbin is a registered trademark of Gyrus ENT LLC.

HydroxylVent is a registered trademark of Gyrus ENT LLC, U.S. Patent No. 5,047,053.

Micron is a registered trademark of Gyrus ENT LLC, U.S. Patent No. 4,744,792.

Hydromer is a trademark of Hydromer, Inc. and is licensed under various US and foreign patents.